

Output Example of the Report Template Files

Output example of the report template files can be shown by clicking the following file name.

(GC/LC Common)

#	File name	Description
1	Method.lsr	Output the method parameters.
2	CalibrationCurveReport_1.lsr	Output the calibration curve graph and the list.
3	CalibrationCurveReport_2.lsr	Output the calibration curve graph and the list in detail.
4	CalibrationCurveReport_3.lsr	Output two blocks of calibration curve graph and the list.
5	PeakReport_1.lsr	Output the chromatogram and the peak report.
6	PeakReport_2.lsr	Output the chromatogram (Landscape) and the peak report.
7	PeakReport_3.lsr	Output the method parameters, the chromatogram and the peak report.
8	PeakReport_4.lsr	Output the chromatogram and the peak report in landscape page.
9	PeakReport_5.lsr	Output the divided chromatogram and the peak report.
10	PeakReport_CalibrationCurve.lsr	Output the chromatogram, the peak table and the calibration curve.
11	QuantitativeResults_1.lsr	Output the all chromatogram and the quantitative results.
12	QuantitativeResults_2.lsr	Output the chromatogram of detector1 and detector2 in a stack and the quantitative results.
13	QuantitativeResults_3.lsr	Output the chromatogram of detector1 and detector2 which are arranged to left and right and the quantitative results.
14	QuantitativeResults_4.lsr	Output the detector1 and the quantitative results, detector2 and the quantitative results.
15	QuantitativeResults_5.lsr	Output the all chromatogram (overlay) and the quantitative results in landscape page.
16	GroupResults_1.lsr	Output the chromatogram and the grouping results.
17	GroupResults_2.lsr	Output the chromatogram, the grouping results and the group calibration curve.
18	ColumnPerformance.lsr	Output the chromatogram and the column performance report.
19	SummaryReport_1.lsr	Output the chromatogram and the peak report for each data.
20	SummaryReport_2.lsr	Output the quantitative summary report for each compound.
21	SummaryReport_3.lsr	Output the concentration summary report, the area summary report and the height summary report.
22	SummaryReport_4.lsr	Output the chromatogram(Landscape), and the quantitative summary report for each compound.
23	SummaryReport_5.lsr	Output the concentration summary report, the area summary report and the height summary report in landscape page.
24	SummaryReport_6.lsr	Output the chromatogram(Overlay), and the quantitative summary report for each compound.
25	BatchTable.lsr	Output the batch table in landscape page.

(LC)

#	File name	Description
1	FractionCollectionReport.lsr	Output the chromatogram and the fraction collection report.
2	PDADataAnalysisResults.lsr	Output the PDA chromatogram, and the spectrum, the library search results and purity of the ID peaks.
3	Contour3DReport.lsr	Output the contour of PDA and the 3D graph.