


# ==== Shimadzu LabSolutions Analysis Report ====

Sample Name : Liquor  
Vial# : 1-4  
Sample Type : Unknown  
Date Acquired : 4/9/2009 2:33:04 AM  
Date Processed : 4/2/2010 4:07:45 PM  
Sample ID : Unknown01  
Injection Volume : 1 uL  
Acquired by : System Administrator  
Modified by : System Administrator

Data Filename : C:\LabSolutions\Sample\GC\Tutorial\_Unk001.gcd  
Method Filename : C:\LabSolutions\Sample\GC\Tutorial\_Method.gcm  
Batch Filename : C:\LabSolutions\Sample\GC\Tutorial Batch.qcb



FID1

ID#	Name	Ret. Time	Area	Height	Conc.	Unit
1	1-Propylalcohol	5.733	20659	8448	85.212	ppm
2	Isobutylalcohol	6.802	45081	17880	143.725	ppm
3	Isoamylalcohol	8.722	90812	36867	277.477	ppm
Total			156551	63195		